

# Holter System



2003 &



2012



Lightweight and compact design digital recorder combined with fast and powerful analysis software, adjustable clinical parameters to meet different practice needs.

- Weighting only 42g
- 3,12 channel digital recording
- Real-time waveform display
- One battery for up to 7 days monitoring
- Removable SD Card
- Pacemaker Detection



**Pulse Healthcare Technology**

House: 11/1, Shahid Minar Road  
Kallyanpur, Dhaka-1216, Bangladesh.

Mobile: +8801708008061

e-mail: [info@pht.com.bd](mailto:info@pht.com.bd)

web: [www.pht.com.bd](http://www.pht.com.bd)

# Holter System



## Holter Software

- Full Disclosure 24/48 Hours or 7 days Recording Time
- Automatic and Manual Editing of Patient Data
- QRS template classification
- Arrhythmia Analysis
- Analysis of ST segment changes and evaluation of myocardial ischemia
- Pacemaker analysis
- Obstructive Sleep Apnea Syndrome analysis
- Atrial Fibrillation analysis
- Review and Print Trends
- Heart Rate Turbulence(HRT) analysis
- Download Data Via Compact Flash Card or USB Cable
- Open communication interfaces of GDT/XML



## Product Specifications

	 2003	 2003 (7 days)	 2012
Recording Channel	3CH	3CH	3, 12 CH
Recording	Full Disclosure, No Data Compression		
Sampling Rate	128Hz, 10,000 s/sec/channel used for pacemaker detection		
A/D Conversion	8 Bit		
Gain	5, 10, 20 mm/mv		
Frequency Response	0.05 Hz to 60 Hz		
Battery	1 AAA alkaline battery, 1.5 Volt		
Operating Duration	24~96 hours	24~168 hours	24~48 hours
Storage	SD Card ( 2G )		
Display	Pacemaker Detection, ECG Display, Lead Quality Check, Low Battery		
Dimensions	2.68"×2.08"×0.63" (68mm×53mm×16mm)		
Weight	1.48 oz (42g)		
Storage	Temperature: -20℃~+60℃ Humidity: 5~95% (non-condensing)		
Operation	Temperature : 0℃~+45℃ Humidity: ≤80% (non-condensing)		



**Pulse Healthcare Technology**  
 House: 11/1, Shahid Minar Road  
 Kallyanpur, Dhaka-1216, Bangladesh.  
 Mobile: +8801708008061  
 e-mail: [info@pht.com.bd](mailto:info@pht.com.bd)  
 web: [www.pht.com.bd](http://www.pht.com.bd)



**ESSE 3** Via Garibaldi 30  
 14022Castelnuovo D.B. (AT)  
 tel +39 011 99 27 706  
 fax +39 011 99 27 506  
 e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
 web : [www.4ci.com](http://www.4ci.com)

