



APEL CO., LTD. JAPAN

**APEL Total Bilirubin Meter
MADE IN JAPAN**

BR-5200P

Hassle-free Intelligent Bilirubin Meter



HIGHLY USER-FRIENDLY

Switch on and Start measuring after 60 sec. warm-up

Very quick measurement

Max. one second measurement time.

No lamp replacement

BR-5200P does not need a lamp replacement thanks to the LED light source.

Compact & light weight

BR-5200P's foot print is small and weights only 1.7kg.

Built-in printer

Compact thermal printer is built in.

Computer connection

BR-5200P has a USB interface for connection with computers.

Free of calibration hassle

Once BR-5200P is calibrated, no more daily calibration is required.



BR-5200P can be connected to a computer via a USB cable.



Measurement data appears on the display.

BR-5200P

Hassle-free Intelligent Bilirubin Meter



**State-of-the-Art
Dual Wavelength
Total Bilirubin Meter
for Neonates**

Ready for measurement:



1 Turn on the switch on the rear panel.



2 Ready for measurement.



3 Insert a capillary tube (containing a sample) into the holder.



4 The value appears very quickly (within one second).

Specifications:

| | |
|--------------------|----------------------------|
| Method | Dual Wavelength |
| Measurement Range | 0-30mg/dL(Total bilirubin) |
| Measuring accuracy | ±5% |
| Sample Volume | 50-60 μ L |
| Display | Large LCD display |
| Sample Container | Hematocrit capillary tube |
| Light source | White high bright LED |
| Detector | Silicon photocell |
| Printer | Thermal serial dot method |
| Interface | USB |
| Power | DC12V 4.5W MAX |
| Dimensions | 262(W) X 214(D) X 109(H)mm |
| Weight | 1.7kg (main body) |

Standard Accessories:

| Description | Quantity |
|-----------------------------------|----------------|
| Capillary tube (Heparinized, Red) | 1 vial(100pcs) |
| Capillary tube (Plain, Blue) | 1 vial(100pcs) |
| Clay sealer (Wax Plate) | 1 pc |
| Standard solution (ABY) | 2 bottles |
| USB cable | 1 pc |
| USB driver software (CD-ROM) | 1 pc |
| APEL AC Adapter set | 1 set |
| Operation manual | 1 copy |

Specifications and appearance are subject to change without notice.

MADE IN JAPAN



APEL CO., LTD. JAPAN

Website

<http://www.apel.co.jp>

Distributor



Pulse Healthcare Technology

House: 11/1, Shahid Minar Road
Kallyanpur, Dhaka-1216, Bangladesh.

Mobile: +8801708008061

e-mail: info@pht.com.bd

web: www.pht.com.bd